# LAMBDA 265

## Users Guide

### Release History

Part Number	Release	Publication Date
09931225	D	September 2016

Any comments about the documentation for this product should be addressed to:

User Assistance  
PerkinElmer  
710 Bridgeport Avenue  
Shelton, CT 06484-4794  
U.S.A.

Or emailed to: info@perkinelmer.com

### Notices

The information contained in this document is subject to change without notice.  
**Except as specifically set forth in its terms and conditions of sale, PerkinElmer makes no warranty of any kind with regard to this document, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose.**  
PerkinElmer shall not be liable for errors contained herein for incidental consequential damages in connection with furnishing, performance or use of this material.

### Copyright Information

This document contains proprietary information that is protected by copyright.  
All rights are reserved. No part of this publication may be reproduced in any form whatsoever or translated into any language without the prior, written permission of PerkinElmer, Inc.

Copyright © 2016 PerkinElmer, Inc.

### Trademarks


Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are protected by law.

PerkinElmer is a registered trademark of PerkinElmer, Inc.

## *Table of Contents*

---

<b>Introduction .....</b>	<b>5</b>
Overview .....	6
About This Manual.....	7
Basic Theory .....	8
PDA UV-Visible Spectrophotometer .....	9
Safety Summary.....	10
Instrument Shipping Kit .....	13
Conventions Used in this Manual.....	14
Notes, cautions and warnings.....	14
Note, attention et avertissement .....	14
Contact Us .....	17
<b>Safety Practices .....</b>	<b>19</b>
Overview .....	20
Precautions.....	21
General Operating Conditions .....	22
Electrical Safety.....	23
EMC Compliance.....	25
Warning Labels .....	26
WEEE Instructions for PerkinElmer Products .....	27
Decontamination.....	28
Cleaning the Instrument .....	28
Fuses .....	29
<b>Lambda 265 Installation.....</b>	<b>31</b>
Preparing the Laboratory.....	32
Installation.....	34
<b>UV Lab Software Installation.....</b>	<b>35</b>
How to Install the UV Lab Software.....	36
<b>Interface Setup .....</b>	<b>41</b>
Setting USB Serial Ports .....	46
<b>Operating Procedure.....</b>	<b>51</b>
Procedure .....	52
<b>Accessories.....</b>	<b>54</b>
Accessories for the Lambda 265.....	55
Single Cell Holder.....	55
Water Jacketed Single Cell Holder.....	56
Variable Pathlength Holder.....	56
Advanced Transmission Holder.....	57
Reflectance Module .....	57
Peltier Controlled Single Cell, With Peltier Controller .....	57
8 Position Multi Cell Holder .....	58
Water Jacketed 8 Position Multi Cell Holder .....	58
Autosampler .....	58
Nano Stick.....	59
Magnetic Stirrer Assembly.....	59
Index .....	60

## *Introduction*

## *Overview*

---

The Lambda 265 instrument runs using UV Lab software.

This manual provides instructions for installing and operating the Lambda 265 UV-Visible Spectrophotometer and UV Lab software. It also introduces and describes how to use basic features of the UV Lab software. For more detailed information on software features including kinetics, bio analysis, color analysis and MCA, refer to the UV Lab software user guide.

## ***About This Manual***

---

This user's guide is divided into following chapters:

### ***Chapter 1 Introduction***

This chapter contains a brief introduction on the instrument, the conventions and warnings used in the manual.

### ***Chapter 2 Safety Practices***

Important safety information is provided in this chapter.

### ***Chapter 3 Lambda 265 Installation***

Information on installing and re-installing your instrument should you ever need to move your system is provided.

### ***Chapter 4 UV Lab Software Installation***

The complete UV Lab software installation procedures for the Lambda 265 is in this chapter.

### ***Chapter 5 Interface Setup***

Interface procedures for the Lambda 265.

### ***Chapter 6 Operating Procedure***

This chapter contains information on operating the Lambda 265.

### ***Chapter 7 Accessories***


This chapter provides information on the accessories available for the Lambda 265.

### ***Index***

## Basic Theory

---


Electromagnetic radiation can pass through most organic and inorganic compounds. Generally, the UV range is from 190 to 380 nm and visible range is from 380-800 nm.

**Figure 1. Electromagnetic Radiation**


When continuous radiation passes through a sample, some is absorbed by the sample. An absorption spectrum can be obtained by monitoring the radiation that penetrates the sample and reaches a detector. For every substance, the absorption rate varies depending on the wavelength of the radiation. Data from absorption spectra can be used for qualitative and quantitative analysis.

A molecule remains in its ground state when it is stable, but it can transition to an excited state when light energy is absorbed. This is termed absorption. (Figure. 2) When the excited molecule returns to the ground state, it emits heat, radiation, fluorescence or phosphorescence. This is termed emission. A molecule whose functional group has a double bond between carbons or between carbon and another atom undergoes transitions in the UV-Visible range. Functional groups that absorb light energy are called Chromophores.

**Figure 2. Energy Absorption**

Once a sample absorbs light energy, the intensity of incoming light  $P_0$  (the energy from the light source) reduces to  $P$  by losing some of its energy; therefore, we define the transmittance of a sample as follows:

**Figure 3. Sample Absorption**

As a result of absorption by a sample, transmittance is defined as;

$$T = P / P_0$$

And Absorbance is defined as:


$$A = -\log T = \log(P_0 / P)$$

Absorbance is related to  $c$  (concentration) and  $b$  (pathlength of light) according to the Lambert-Beer law and the equation;

$$A = \epsilon bc$$

Where,  $\epsilon$ , the molar absorption coefficient of a particular molecule, varies with the wavelength.

Although Beer's law has several limitations, it is widely applied to quantitative analysis.

**Figure 4. Transmittance and Absorbance**

## **PDA UV-Visible Spectrophotometer**

The PDA Spectrophotometer has a multi-channel detector controlled by a microprocessor which collects spectral data for several wavelengths simultaneously. In contrast, a PMT Spectrophotometer has a single channel detector which collects data for each specific wavelength in turn. Therefore the PDA Spectrophotometer is much faster and gives more reproducible results than the PMT Spectrophotometer which depends on a rotating grating or prism operated by a stepper motor to obtain individual wavelengths.

One feature of a PDA is an integrating photo-detector, which integrates charge depending on the light intensity. The advantage of the integrating function is named Felgette's S/N advantage or Multi-channel advantage.

## Safety Summary


---

The following safety symbols are used on this product.

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"OFF" (power)

To indicate disconnection from the mains, at least for mains switches or their positions, and all those case where safety is involved.

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The user cannot repair this equipment.

DO NOT attempt to open the case or to disassemble internal parts.

Only your PerkinElmer service representative with knowledge of electrical fire and shock hazard should service this equipment.

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This equipment should be used in the correct position. If the equipment is turned on its side, etc., it will be unstable and may be damaged if it falls over as a result of receiving a slight mechanical shock.

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## ***Instrument Shipping Kit***

---

The following table shows what is included in the Lambda 265 Spectrophotometer shipping kit.

<b>Item</b>	<b>Name</b>	<b>Qty.</b>
Main Instrument	Lambda 265 Spectrophotometer	1
Accessories	Power Cord (US/EU/UK)	1/1/1
	Interface Cable (USB)	1
	Fuse AC 250V (T2.0AL)	2
	Wrench, 3mm	1
	Document Pack CD	1
	UV Lab Software CD	1
	Single Cell Holder	1

---

## ***Conventions Used in this Manual***

---

Normal text is used to provide information and instructions.

**Bold** text refers to text that is displayed on the screen.

UPPERCASE text, for example ENTER or ALT, refers to keys on the PC keyboard. '+' is used to show that you have to press two keys at the same time, for example, ALT+F.

All eight digit numbers are PerkinElmer part numbers unless stated otherwise.

### ***Notes, cautions and warnings***

Three terms, in the following standard formats, are also used to highlight special circumstances and warnings.

**NOTE:** *A note indicates additional, significant information that is provided with some procedures.*

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**WARNING**

We use the term **WARNING** to inform you about situations that could result in **personal injury** to yourself or other persons. Details about these circumstances are in a box like this one.

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## ***Contact Us***

---

Supplies, replacement parts, and accessories can be ordered directly from PerkinElmer, using the part numbers quoted in the guides provided with the instrument.

See our website:


<http://perkinelmer.com>

PerkinElmer's catalog service offers a full selection of high-quality supplies.

To place an order for supplies and many replacement parts, request a free catalog, or ask for information:

If you are located within the U.S., call toll free 1-800-762-4000, 8 a.m. to 8 p.m. EST. Your order will be shipped promptly, usually within 24 hours.

If you are located outside of the U.S., call your local PerkinElmer sales or service office.

## *Safety Practices*

## ***Overview***

---

This chapter describes the general safety practices and precautions that must be observed when operating the Lambda 265.

This advice is intended to supplement, not supersede, the normal safety codes in the user's country. It is also a supplement to the PerkinElmer standard Safety and Health Policy. The information provided does not cover every safety procedure that should be practiced. Ultimately, maintenance of a safe laboratory environment is the responsibility of the analyst and the analyst's organization.

Please consult all manuals supplied with the Lambda 265 and accessories before you start working with the instrument. Carefully read the safety information in this chapter and in the other manuals supplied. When setting up the instrument or performing analyses or maintenance procedures, strictly follow the instructions provided.

## *Precautions*

---

**WARNING**

*Be sure that all instrument operators read and understand the precautions listed below. It is advisable to post a copy of the precautions near or on the instrument shelf.*

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The following precautions must be observed when using the Lambda 265:

- Be sure that the voltage of the Lambda 265 corresponds to the voltage used in your laboratory.

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## General Operating Conditions

---

The Lambda 265 has been designed and tested in accordance with PerkinElmer specifications and in accordance with the safety requirements of the International Electrotechnical Commission (IEC). The Lambda 265 conforms to IEC61010-1 (Safety Requirements for electrical equipment for measurement, control and laboratory use) as it applies to IEC Class 1 (earthed) appliances and therefore meets the requirements of EC directive 2006/95/EC.

Only use the Lambda 265 indoors and under the following conditions:

Temperature            15 °C to 35 °C (59-95 °F)

Relative Humidity      Less than 80% relative non-condensing

Indoor only            Altitude up to 2000m

If possible, avoid any adjustment, maintenance and repair of the opened, operating instrument. If any adjustment, maintenance and repair of the opened instrument is necessary, this must be done by a skilled person who is aware of the hazard involved.

Whenever it is likely that the Lambda 265 is unsafe, make it inoperative. The Lambda 265 may be unsafe if it:

- shows visible damage
- fails to perform the intended measurement
- has been subjected to prolonged storage in unfavourable conditions
- has been subjected to severe transport stresses.

**WARNING**

*If the equipment is used in a manner not specified herein, the protection provided by the equipment may be impaired.*

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## Electrical Safety

---

The instrument has been designed to protect the operator from potential electrical hazards. This section describes some recommended electrical safety practices.

**WARNING**

*Lethal voltages are present at certain areas within the instrument. Installation and internal maintenance of the instrument should only be performed by a PerkinElmer service engineer or similarly authorized and trained person. When the instrument is connected to line power, opening the instrument covers is likely to expose live parts. Even when the power switch is off, high voltages can still be present. Capacitors inside the instrument may still be charged even if the instrument has been disconnected from all voltage sources.*

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The instrument must be correctly connected to a suitable electrical supply. The supply must have a correctly installed protective conductor (earth ground) and must be installed or checked by a qualified electrician before connecting the instrument.

**WARNING**

*Any interruption of the protective conductor (earth ground) inside or outside the instrument or disconnection of the protective conductor terminal is likely to make the instrument dangerous. Intentional interruption is prohibited.*

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When working with the instrument:

- Connect the instrument to a correctly installed line power outlet that has a protective conductor connection (earth ground).
- Do not operate the instrument with any covers or internal parts removed.
- Do not attempt to make internal adjustments or replacements except as directed in the manuals.

- Disconnect the instrument from all voltage sources before opening it for any adjustment, replacement, maintenance, or repair. If afterwards, the opened instrument must be operated for further adjustment, maintenance, or repair, this must only be done by a PerkinElmer Service engineer.
- Whenever it is possible that the instrument is no longer electrically safe for use, make the instrument inoperative and secure it against any unauthorized or unintentional operation. The electrical safety of the instrument is likely to be impaired if, for example, the instrument shows visible damage; has been subjected to prolonged storage under unfavorable conditions; or has been subjected to severe stress during transportation.

## ***EMC Compliance***

---

### ***EC directive***

The Lambda 265 has been designed and tested to meet the requirements of the EC directive 2004/108/EC. The Lambda 265 complies with the EMC standard EN61326, (EMC standard for electrical equipment for measurement, control and laboratory use).

### ***FCC rules and regulations***

This product is classified as a digital device used exclusively as industrial, commercial, or medical test equipment. It is exempt from the technical standards specified in Part 15 of the FCC Rules and Regulations based on Section 15.103 (c).

### **Installation Category**

This product will operate safely in electrical power environments that are rated as Installation Category II in EN/IEC 61010-1.

### **Pollution Degree**

This product will operate safely in environments that contain nonconductive foreign matter up to Pollution Degree 2 in EN/IEC 61010-1.

## Warning Labels

**WARNING**

*When this label is attached to an instrument it means refer to the manual.*

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See the following figure for the location of warning labels on the Lambda 265.

Warning labels on the back of the instrument:


**Warning  
FUSE**

**Avertissement  
FUSE**

## ***WEEE Instructions for PerkinElmer Products***

---

A label with a crossed-out wheeled bin symbol and a rectangular bar indicates that the product is covered by the Waste Electrical and Electronic Equipment (WEEE) Directive and is not to be disposed of as unsorted municipal waste. Any products marked with this symbol must be collected separately, according to the regulatory guidelines in your area.

The objectives of this program are to preserve, protect and improve the quality of the environment, protect human health, and utilize natural resources prudently and rationally. Specific treatment of WEEE is indispensable in order to avoid the dispersion of pollutants into the recycled material or waste stream. Such treatment is the most effective means of protecting the customer's environment.

Requirements for waste collection reuse, recycling, and recovery programs vary by regulatory authority at your location. Contact your local responsible body (e.g., your laboratory manager) or authorized representative for information regarding applicable disposal regulations. Contact PerkinElmer at the web site listed below for information specific to PerkinElmer products.

**Web address:** <http://www.perkinelmer.com/pages/010/onesource/environmental-health-and-safety/environmental-directives-compliance.xhtml> For Customer Care telephone numbers select "Contact us" on the web page.

Products from other manufacturers may also form a part of your PerkinElmer system. These other producers are directly responsible for the collection and processing of their own waste products under the terms of the WEEE Directive. Please contact these producers directly before discarding any of their products.

Consult the PerkinElmer web site (above) for producer names and web addresses.

## ***Decontamination and Cleaning***

---

Before using any cleaning or decontamination methods except those specified by PerkinElmer, users should check with PerkinElmer that the proposed method will not damage the equipment.

### ***Decontamination***

Customers wishing to return instrumentation and/or associated materials to PerkinElmer for repair, maintenance, warranty or trade-in purposes are advised that all returned goods must be certified as clean and free from contamination.

The customer's responsible body is required to follow the "Equipment Decontamination Procedure" and complete the "Certificate of Decontamination". These documents are available on the PerkinElmer public website:

**Procedure:**

[http://las.perkinelmer.com/Content/technicalinfo/dts\\_instrumentdeconprocedure.pdf](http://las.perkinelmer.com/Content/technicalinfo/dts_instrumentdeconprocedure.pdf)

**Certificate form:**

[http://las.perkinelmer.com/Content/technicalinfo/dts\\_perkinelmercertificationofdecontaminationform.pdf](http://las.perkinelmer.com/Content/technicalinfo/dts_perkinelmercertificationofdecontaminationform.pdf)

If you do not have access to the internet and are located in the U.S., call toll free at 1-800-762-4000 or (+1) 203-925-4602, 8:30 a.m. - 7 p.m. EST and speak to Customer Support.

In Canada, call toll free 800-561-4646 and speak to Customer Support.

If you are located outside of the United States or Canada, please call your local PerkinElmer sales office for more information.

### ***Cleaning the Instrument***

Exterior surfaces may be cleaned with a soft cloth, dampened with a mild detergent and water solution. Do **not** use abrasive cleaners or solvents.

## Fuses

*There is risk of receiving a fatal electric shock if the fuses are replaced with the power cord connected.*

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
The fuse is located at the rear of the instrument.

1. Turn off and unplug the instrument.
2. Carefully open the compartment latch where the fuse is located.

Location of latch on the  
fuse compartment door

3. Disconnect the fuse.
4. Replace with a new T2.0AL fuse. One spare is contained in the power module.
5. Close the compartment door.
6. Plug in the instrument and turn on.

## Lambda 265 Installation

## ***Preparing the Laboratory***

---

The compact size of the Lambda 265 enables you to install the spectrophotometer in a limited space.

### ***Make Room for installation***

The Lambda 265 System needs sufficient space for a computer, printer and selected accessories such as the Auto Sampler. The system dimensions are 340 mm (W) × 320 mm (D) × 115 mm (H) [13.4 in. (W) x 12.6 in. (D) x 4.5 in. (H)].

**WARNING**

*In view of the electronic components in the Lambda 265 System, It should be positioned well away from any instruments emitting an electromagnetic field and any environment that may cause vibration.*

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### ***Power***

A grounded AC power supply must be used. The Lambda 265 accepts 100 to 240 VAC at 50 Hz or 60 Hz.

**WARNING**

*In order to protect the system from electrical noise and current surges, please connect the power cord to a GROUNDED OUTLET.*

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### ***Environmental Requirements***

This system is for indoor use only.

<b>Requirement</b>	<b>Description</b>
Indoor use	Altitude up to 2000m (6,500 ft.)
	Pollution degree 2
Dimensions	340 mm(W) × 320 mm(D) × 115 mm(H)
Weight	6.6 kg (14.6 lbs)
Electrical supply	AC 100-240V, 50/60Hz, 20W
Operating Temperature	15 - 35°C (59 - 95°F)
Operating humidity	Less than 80% relative, non-condensing

### ***Minimum Computer Requirements .***

<b>Requirement</b>	<b>Description</b>
Hardware	PC (Intel® Core 1.5GHz or faster processor)  At least 1 GB RAM  50 GB with 1 GB free HDD  Color Monitor Display (1024 × 768)  CD-ROM Drive  Mouse & Keyboard  USB port for the data acquisition
Operating System	Microsoft® Windows 7
Output device	Microsoft® Windows compatible Printer

## Installation

---

1. Carefully remove the instrument and all of its components from the shipping container.
2. Verify that all of the components received match the packing list.
3. Verify that all system components required for installation are present, this, must include:
  - Lambda 265 PDA UV-Visible Spectrophotometer
  - Interface cable (USB)
  - Power Cord
  - UV Lab Software CD
  - Computer systems
4. Place the instrument in a location that is compatible with the required environmental conditions for operation. Connect the power cord to the plug; the power plug must be grounded.


Figure 4 System in location

**WARNING**

*The electric power could be dangerous to the body. Before connecting the power cord, verify that the main power switch is off.*

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5. Setup the up the computer using the manufacture's supplied documentation.

In the following chapters you will:


- Install UV Lab software.
- Install the USB driver for Interface with the computer.
- Execute UV Lab software.
- Installation is complete and the system is ready for operation.

## *UV Lab Software Installation*


## How to Install the UV Lab Software

---


1. Insert the UV Lab installation CD into the CD-ROM drive. InstallShield® Wizard starts automatically.
2. InstallShield® Wizard does not start automatically, select **Computer** → **CD Drive**. Double click **setup.exe**.


3. UV Lab InstallShield Wizard starts. Select **Next**.


4. Enter customer information, if necessary, and select **Next**.


5. The software will be installed in C:\ UV Lab as a destination folder. To install the software to a different folder, select **Change** and select the desired folder.


6. Verify the current settings and select **Install**.


7. The following window shows the status of the UV Lab installation.

8. After the setup is completed successfully, select **Finish**. The UV Lab software icon will appear on the Desktop.


## Interface Setup


## Installing the USB Drive

---


1. Turn on the computer. And then insert the enclosed CD into the CD-ROM drive.
2. Connect the USB between the computer and instrument.
3. Click **Computer** → **CD Drive** → **Open**.


4. Click **USB Drivers** folder.


5. Select **LAMBDA 265** folder.


6. Select **Windows XP\_7\_8, Vista, Server 2003 to 2008 win2000 R2** folder.


7. Double Click **USBCom-CDM\_20828**.


8. Select **Extract**.

9. Select **Next**.


10. The following dialog box will appear. After installation is completed successfully, select **Finish**.


11. Change the USB serial Ports after finishing driver setup. See the following chapter.

## Setting USB Serial Ports


1. Select Computer → Properties.


2. Select on Device Manager.


3. Select the **Ports (COM & LPT)** to expand the listing. These are the devices currently connected to the COM ports. **USB Serial Port (COMx)** is visible when the driver installation is completed successfully.


4. Double click **USB Serial Port (COMx)** of the Ports (COM & LPT) section.


5. Select on the **Port Settings** tab and select on the **Advanced...** button.


6. Change the parameter values as shown below.


7. Select **OK** after checking the changed parameter values.

8. Launch the **UV Lab** software.
9. If the Lambda 265 fails to communicate with the PC, change the COM Port Number as following procedure.
10. Open **Advanced Setting for COMx** window following procedure on page 469.
11. Select on the **COM Port number** list to expand it and change the COM port number another one which is not in use from COM 1 to COM 10.

12. Make sure that the changed COM Port Number is applied, and select **OK**.


13. Restart the computer after finishing driver setting.

## *Operating Procedure*


## Procedure

---

1. Turn on the power switch.
2. Start **UV Lab** software.
3. Enter a Title, Sample Name and choose the **Experiment Type**. Select **OK**.


4. Set method parameters. Select **OK**.


5. Prepare blank and sample solutions.
6. Fill a 10 mm pathlength cell to ~80% of its volume with the blank solution and place it into the cell holder. Select or select **Run Blank** from the **Measure** menu.

7. The blank spectrum is displayed.
8. Remove the cell from the cell holder.
9. Fill a 10 mm pathlength cell to ~80% of its volume with the sample solution and place it into the cell holder. Click  Sample or select Run Sample in the Measure menu.
10. The sample spectrum is displayed.

11. Print and Save the spectral data.

## Accessories

## *Accessories for the Lambda 265*

---

Every day you can count on PerkinElmer to provide you with solutions that deliver reliable performance, control operating costs and maximize operational time. Our complete portfolio of consumables, parts, supplies, training and service helps you meet both routine and demanding measurement challenges. We invest heavily in testing and validating our products to ensure you receive guaranteed compatibility and performance-on-time, every time, for every instrument in your laboratory.

Supplies, replacement parts, and accessories can be ordered directly from PerkinElmer, using the part numbers quoted in the guides provided with the instrument.

See our website:

<http://perkinelmer.com>

PerkinElmer's catalog service offers a full selection of high-quality supplies.

To place an order for supplies and many replacement parts, request a free catalog, or ask for information:

If you are located within the U.S., call toll free 1-800-762-4002, 8 a.m. to 8 p.m. EST. Your order will be shipped promptly, usually within 24 hours.

If you are located outside of the U.S., call your local PerkinElmer sales or service office.

### *Single Cell Holder*

Designed for a standard cell with a 10 mm pathlength. Optional stirring is available with a separately purchased stir controller.

## ***Water Jacketed Single Cell Holder***

Use this accessory to control the temperature of a sample in a standard 10 mm rectangular cuvette. Optional stirring is available with a separately purchased stir controller. This accessory accommodates a temperature by circulator.

## ***Variable Pathlength Holder***

Use this accessory to hold rectangular cells with 5, 10, 20, 40, 50 mm pathlengths. A moveable block secures cells in position.

## *Advanced Transmission Holder*

Use this accessory to measure the transmittance spectrum of various solid samples including films, glass, filters and optical components.

## *Reflectance Module*

Use this accessory to test coatings and semi-conductor wafers; especially effective for measuring film thickness.

## *Peltier Controlled Single Cell, With Peltier Controller*

Control the temperature of a single 10 mm pathlength rectangular cell between -5°C and 100°C. This accessory includes a temperature controller, Peltier accessory and all accessories necessary for its operation and connection to the instrument.

## **8 Position Multi Cell Holder**

Use this accessory to measure eight samples automatically. This accessory is controlled by UV Lab software.

## **Water Jacketed 8 Position Multi Cell Holder**

Use this accessory to measure eight samples automatically. This accessory accommodates a temperature by circulator. This accessory is controlled by UV Lab software.

## **Autosampler**

Automate the measurement of up to 360 samples at a time with an autosampler accessory. Two and four rack configurations are available. This accessory system, including the automated sippier, is fully controlled by UV Lab software.

## ***Nano Stick***

Use this accessory to measure micro-volume sample such as DNA, RNA and Protein.

## ***Magnetic Stirrer Assembly***

Use this accessory to mix sample, this accessory is controlled by Control Unit. Control the speed of magnetic stirrer between 120 and 1200 rpm.

## ***Index***

---

<b>A</b>	SB drive ..... 42
Accessories ..... 58	
<b>B</b>	
Basic Theory ..... 8	
<b>C</b>	
Cleaning ..... 28	
Cleaning the Instrument ..... 28	
Computer Requirements ..... 33	
Conventions	
Note, attention et avertissement ..... 14	
Notes, cautions and warnings ..... 14	
text ..... 14	
<b>D</b>	
Decontamination ..... 28	
<b>E</b>	
Electrical Safety ..... 23	
Environmental Requirements ..... 33	
<b>I</b>	
Installation	
Lambda 265 ..... 34	
<b>L</b>	
Laboratory Preparation ..... 32	
<b>O</b>	
Operating Conditions ..... 22	
Operating Procedure ..... 54	
<b>P</b>	
Precautions ..... 21	
<b>S</b>	
Safety	
Summary ..... 10	
Safety Information ..... 21	
electrical safety ..... 23	
Shipping Kit ..... 13	
<b>U</b>	
USB drive	
installation ..... 42	
USB Serial Ports	
setup ..... 47	