LAMBDA 365

Users Guide

Release History

| Part Number | Release | Publication Date |
|-------------|---------|------------------|
| 09931249 | E | September 2016 |

Any comments about the documentation for this product should be addressed to:

User Assistance
PerkinElmer
710 Bridgeport Avenue
Shelton, CT 06484-4794
U.S.A.

Or emailed to: info@perkinelmer.com

Notices

The information contained in this document is subject to change without notice.

Except as specifically set forth in its terms and conditions of sale, PerkinElmer makes no warranty of any kind with regard to this document, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose.

PerkinElmer shall not be liable for errors contained herein for incidental consequential damages in connection with furnishing, performance or use of this material.

Copyright Information

This document contains proprietary information that is protected by copyright.

All rights are reserved. No part of this publication may be reproduced in any form whatsoever or translated into any language without the prior, written permission of PerkinElmer, Inc.

Copyright © 2015 PerkinElmer, Inc.


Trademarks

Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are protected by law.

PerkinElmer is a registered trademark of PerkinElmer, Inc.

Table of Contents

| | |
|--|-----------|
| Table of Contents | 3 |
| Introduction | 5 |
| Overview | 6 |
| About This Manual | 7 |
| Basic Theory | 8 |
| Double Beam UV-Visible Spectrophotometer | 10 |
| Safety Summary | 11 |
| Instrument Shipping Kit | 14 |
| Conventions Used in this Manual | 15 |
| Notes, cautions and warnings | 15 |
| Note, attention et avertissement | 15 |
| Contact Us | 18 |
| Safety Practices | 19 |
| Overview | 20 |
| Precautions | 21 |
| General Operating Conditions | 22 |
| Electrical Safety | 23 |
| EMC Compliance | 25 |
| Warning Labels | 26 |
| WEEE Instructions for PerkinElmer Products | 27 |
| Decontamination | 28 |
| Cleaning the Instrument | 28 |
| Fuses | 29 |
| Lambda 365 Installation | 31 |
| Preparing the Laboratory | 32 |
| Installation | 34 |
| UV Express Software Installation | 35 |
| How to Install the UV Express Software | 36 |
| Interface Setup | 41 |
| Setting USB Serial Ports | 47 |
| Procedure | 48 |
| Operating Procedure | 53 |
| Procedure | 54 |
| Accessories | 55 |
| Accessories for the Lambda 365 | 56 |
| Single Cell Holder | 56 |
| Water Jacketed Single Cell Holder | 57 |
| 8 Position Multicell Holder / Water Jacketed 8-Position Multicell Holder | 57 |
| Micro Cell Holder | 57 |
| Variable Pathlength Cuvette Holder | 58 |
| Advanced Transmission Holder | 58 |
| Film Holder | 58 |
| Variable Angle Transmission Holder | 59 |
| 50 mm Transmission / Reflectance Sphere | 59 |
| Fixed Angle Reflectance Holder | 59 |
| Single Cell Peltier Holder, Controller | 60 |
| 6-Position Peltier Controlled Cell Changer, Controller | 60 |
| Nano Stick | 61 |
| Autosampler | 61 |
| Lamp Replacement | 63 |
| D2 Lamp Replacement | 64 |
| Tungsten Lamp Module Replacement | 66 |
| Index | 70 |

Introduction

Overview

This manual provides instructions for installing and operating the Lambda 365 Double Beam UV-Visible Spectrophotometer and UV Express software. For more detailed information on software features including Quantification, kinetics and bio software, refer to the UV Express software user guide.

About This Manual

This user's guide is divided into following chapters:

Chapter 1 Introduction

This chapter contains a brief introduction on the instrument, the conventions and warnings used in the manual.

Chapter 2 Safety Practices

Important safety information is provided in this chapter.

Chapter 3 Lambda 365 Installation

Information on installing and re-installing your instrument should you ever need to move your system is provided.

Chapter 4 UV Express Software Installation

The complete UV Express software installation procedures for the Lambda 365 is in this chapter.

Chapter 5 Interface Setup

Interface procedures for the Lambda 365.

Chapter 6 Operating Procedure

This chapter contains information on operating the Lambda 365.

Chapter 7 Accessories

This chapter provides information on the accessories available for the Lambda 365.

Index

Basic Theory

Electromagnetic radiation can pass through most organic and inorganic compounds. Generally, the UV range is from 190 to 380 nm and visible range is from 380-800 nm.


Figure 1. Electromagnetic Radiation

When continuous radiation passes through a sample, some is absorbed by the sample. An absorption spectrum can be obtained by monitoring the radiation that penetrates the sample and reaches a detector. For every substance, the absorption rate varies depending on the wavelength of the radiation. Data from absorption spectra can be used for qualitative and quantitative analysis.

A molecule remains in its ground state when it is stable, but it can transition to an excited state when light energy is absorbed. This is termed absorption. (Figure. 2) When the excited molecule returns to the ground state, it emits heat, radiation, fluorescence or phosphorescence. This is termed emission. A molecule whose functional group has a double bond between carbons or between carbon and another atom undergoes transitions in the UV-Visible range. Functional groups that absorb light energy are called Chromophores.


Figure 2. Energy Absorption

Once a sample absorbs light energy, the intensity of incoming light P_0 (the energy from the light source) reduces to P by losing some of its energy;

therefore, we define the transmittance of a sample as follows:


Figure 3. Sample Absorption

As a result of absorption by a sample, transmittance is defined as;

$$T = P / P_0$$

And Absorbance is defined as:

$A = -\log T = \log(P_0 / P)$ Absorbance is related to c (concentration) and b (pathlength of light) according to the Lambert-Beer law and the equation;

$$A = \epsilon bc$$

Where, ϵ , the molar absorption coefficient of a particular molecule, varies with the wavelength.

Although Beer's law has several limitations, it is widely applied to quantitative analysis.



Figure 4. Transmittance and Absorbance

Double Beam UV-Visible Spectrophotometer

Lambda 365 double-beam UV-Visible spectrophotometer, the light is split into two beams before it reaches the sample. One beam is used as the reference; the other beam passes through the sample. The reference beam intensity is taken as 100% Transmission (or 0 Absorbance), and the measurement displayed is the ratio of the two beam intensities. Lambda 365 double-beam UV-Visible spectrophotometer has two detectors (photodiodes), and the sample and reference beam are measured at the same time.


Safety Summary

The following safety symbols are used on this product.

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"OFF" (power)

To indicate disconnection from the mains, at least for mains switches or their positions, and all those case where safety is involves.

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The user cannot repair this equipment.

DO NOT attempt to open the case or to disassemble internal parts.

Only PerkinElmer service representative with knowledge of electrical fire and shock hazard should service this equipment.

Repair

There are high-voltage parts in this equipment presenting a risk of severe injury or fatal electric shock to untrained personnel.

In addition, there is a risk of damage to precision parts.

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This equipment should be used in the correct position. If the equipment is turned on its side, etc., it will be unstable and may be damaged if it falls over as a result of receiving a slight mechanical shock.

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Instrument Shipping Kit

The following table shows what is included in the Lambda 365 Spectrophotometer shipping kit.

| Item | Name | Qty. |
|-----------------|---|-------------|
| Main Instrument | Lambda 365 Double Beam UV-Vis Spectrophotometer | 1 |
| Accessories | Power Cord (US/EU/UK) | 3 |
| | Interface Cable (USB) | 1 |
| | Lambda 365 Document Pack CD | 1 |
| | UV Express Software CD | 1 |
| | Single Cell Holder | 1 |
| | Lamp Change Tool | 1 |
| | Fuse AC 250 V (T3.15AL) | 2 |
| | Wrench (M3) | 1 |

Conventions Used in this Manual

Normal text is used to provide information and instructions.

Bold text refers to text that is displayed on the screen.

UPPERCASE text, for example ENTER or ALT, refers to keys on the PC keyboard. '+' is used to show that you have to press two keys at the same time, for example, ALT+F.

All eight digit numbers are PerkinElmer part numbers unless stated otherwise.

Notes, cautions and warnings

Three terms, in the following standard formats, are also used to highlight special circumstances and warnings.

NOTE: *A note indicates additional, significant information that is provided with some procedures.*

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**WARNING**

We use the term **WARNING** to inform you about situations that could result in **personal injury** to yourself or other persons. Details about these circumstances are in a box like this one.

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Contact Us

Supplies, replacement parts, and accessories can be ordered directly from PerkinElmer, using the part numbers quoted in the guides provided with the instrument.

See our website:

<http://perkinelmer.com>

PerkinElmer's catalog service offers a full selection of high-quality supplies.

To place an order for supplies and many replacement parts, request a free catalog, or ask for information:

If you are located within the U.S., call toll free 1-800-762-4002, 8 a.m. to 8 p.m. EST. Your order will be shipped promptly, usually within 24 hours.

If you are located outside of the U.S., call your local PerkinElmer sales or service office.

Safety Practices

Overview

This chapter describes the general safety practices and precautions that must be observed when operating the Lambda 365.

This advice is intended to supplement, not supersede, the normal safety codes in the user's country. It is also a supplement to the PerkinElmer standard Safety and Health Policy. The information provided does not cover every safety procedure that should be practiced. Ultimately, maintenance of a safe laboratory environment is the responsibility of the analyst and the analyst's organization.

Please consult all manuals supplied with the Lambda 365 and accessories before you start working with the instrument. Carefully read the safety information in this chapter and in the other manuals supplied. When setting up the instrument or performing analyses or maintenance procedures, strictly follow the instructions provided.

Precautions

WARNING

Be sure that all instrument operators read and understand the precautions listed below. It is advisable to post a copy of the precautions near or on the instrument shelf.

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The following precautions must be observed when using the Lambda 365:

- Be sure that the voltage of the Lambda 365 corresponds to the voltage used in your laboratory.

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General Operating Conditions

The Lambda 365 has been designed and tested in accordance with PerkinElmer specifications and in accordance with the safety requirements of the International Electrotechnical Commission (IEC). The Lambda 365 conforms to IEC61010-1 (Safety Requirements for electrical equipment for measurement, control and laboratory use) as it applies to IEC Class 1 (earthed) appliances and therefore meets the requirements of EC directive 2006/95/EC.

Only use the Lambda 365 indoors and under the following conditions:

| | |
|-------------------|---------------------------------------|
| Temperature | 15 °C to 35 °C |
| Relative Humidity | Less than 80% relative non-condensing |
| Indoor only | Altitude up to 2000 m |

If possible, avoid any adjustment, maintenance and repair of the opened, operating instrument. If any adjustment, maintenance and repair of the opened instrument is necessary, this must be done by a skilled person who is aware of the hazard involved.

Whenever it is likely that the Lambda 365 is unsafe, make it inoperative. The Lambda 365 may be unsafe if it:

- shows visible damage
- fails to perform the intended measurement
- has been subjected to prolonged storage in unfavourable conditions
- has been subjected to severe transport stresses.

WARNING

If the equipment is used in a manner not specified herein, the protection provided by the equipment may be impaired.

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Electrical Safety

The instrument has been designed to protect the operator from potential electrical hazards. This section describes some recommended electrical safety practices.

WARNING

Lethal voltages are present at certain areas within the instrument. Installation and internal maintenance of the instrument should only be performed by a PerkinElmer service engineer or similarly authorized and trained person. When the instrument is connected to line power, opening the instrument covers is likely to expose live parts. Even when the power switch is off, high voltages can still be present. Capacitors inside the instrument may still be charged even if the instrument has been disconnected from all voltage sources.

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The instrument must be correctly connected to a suitable electrical supply. The supply must have a correctly installed protective conductor (earth ground) and must be installed or checked by a qualified electrician before connecting the instrument.

WARNING

Any interruption of the protective conductor (earth ground) inside or outside the instrument or disconnection of the protective conductor terminal is likely to make the instrument dangerous. Intentional interruption is prohibited.

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When working with the instrument:

- Connect the instrument to a correctly installed line power outlet that has a protective conductor connection (earth ground).
- Do **not** operate the instrument with any covers or internal parts removed.

24 . Lambda 365 Users Guide

- Do **not** attempt to make internal adjustments or replacements except as directed in the manuals.
- Disconnect the instrument from all voltage sources before opening it for any adjustment, replacement, maintenance, or repair. If afterwards, the opened instrument must be operated for further adjustment, maintenance, or repair, this must only be done by a PerkinElmer Service engineer.
- Whenever it is possible that the instrument is no longer electrically safe for use, make the instrument inoperative and secure it against any unauthorized or unintentional operation. The electrical safety of the instrument is likely to be impaired if, for example, the instrument shows visible damage; has been subjected to prolonged storage under unfavorable conditions; or has been subjected to severe stress during transportation.

EMC Compliance

EC directive

The Lambda 365 has been designed and tested to meet the requirements of the EC directive 2006/95/EC. The Lambda 365 complies with the EMC standard EN61326-1, (EMC standard for electrical equipment for measurement, control and laboratory use) and EN61000-3.

FCC rules and regulations

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Installation Category

This product will operate safely in electrical power environments that are rated as Installation Category II in EN/IEC 61010-1.

Pollution Degree


This product will operate safely in environments that contain nonconductive foreign matter up to Pollution Degree 2 in EN/IEC 61010-1.

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See the following figure for the location of warning labels on the Lambda 365.

Warning labels on the back of the instrument:

Warning
For continued protection against electrical shock, replace with same type and ratings of fuse.

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Warning
FUSE

Avertissement
FUSE

Warning
To prevent electrical shock, do not open cover.
Refer all servicing to qualified personnel.

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WEEE Instructions for PerkinElmer Products

or

A label with a crossed-out wheeled bin symbol and a rectangular bar indicates that the product is covered by the Waste Electrical and Electronic Equipment (WEEE) Directive and is not to be disposed of as unsorted municipal waste. Any products marked with this symbol must be collected separately, according to the regulatory guidelines in your area.

The objectives of this program are to preserve, protect and improve the quality of the environment, protect human health, and utilize natural resources prudently and rationally. Specific treatment of WEEE is indispensable in order to avoid the dispersion of pollutants into the recycled material or waste stream. Such treatment is the most effective means of protecting the customer's environment.

Requirements for waste collection reuse, recycling, and recovery programs vary by regulatory authority at your location. Contact your local responsible body (e.g., your laboratory manager) or authorized representative for information regarding applicable disposal regulations. Contact PerkinElmer at the web site listed below for information specific to PerkinElmer products.

Web address: <http://www.perkinelmer.com/pages/010/onesource/environmental-health-and-safety/environmental-directives-compliance.xhtml>

For Customer Care telephone numbers select "Contact us" on the web page.

Products from other manufacturers may also form a part of your PerkinElmer system. These other producers are directly responsible for the collection and processing of their own waste products under the terms of the WEEE Directive. Please contact these producers directly before discarding any of their products.

Consult the PerkinElmer web site (above) for producer names and web addresses.

Decontamination and Cleaning

Before using any cleaning or decontamination methods except those specified by PerkinElmer, users should check with PerkinElmer that the proposed method will not damage the equipment.

Decontamination

Customers wishing to return instrumentation and/or associated materials to PerkinElmer for repair, maintenance, warranty or trade-in purposes are advised that all returned goods must be certified as clean and free from contamination.

The customer's responsible body is required to follow the "Equipment Decontamination Procedure" and complete the "Certificate of Decontamination". These documents are available on the PerkinElmer public website:

Procedure:

http://las.perkinelmer.com/Content/technicalinfo/dts_instrumentdeconprocedure.pdf

Certificate form:

http://las.perkinelmer.com/Content/technicalinfo/dts_perkinelmercertificationofdecontaminationform.pdf

If you do not have access to the internet and are located in the U.S., call toll free at 1-800-762-4000 or (+1) 203-925-4602, 8:30 a.m. - 7 p.m. EST and speak to Customer Support.

In Canada, call toll free 800-561-4646 and speak to Customer Support.

If you are located outside of the United States or Canada, please call your local PerkinElmer sales office for more information.

Cleaning the Instrument

Exterior surfaces may be cleaned with a soft cloth, dampened with a mild detergent and water solution. Do **not** use abrasive cleaners or solvents.

Fuses

WARNING

There is risk of receiving a fatal electric shock if the fuses are replaced with the power cord connected.

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The fuse is located at the rear of the instrument.

1. Turn off and unplug the instrument.
2. Carefully open the compartment latch where the fuse is located.

Location of latch
on the fuse
compartment
door

3. Disconnect the fuse.
4. Replace with two new T3.15L fuses, one for each phase.
5. Close the compartment door.
6. Plug in the instrument and turn on.

Lambda 365 Installation

Preparing the Laboratory

The compact size of the Lambda 365 enables you to install the spectrophotometer in a limited space.

Make Room for installation

The Lambda 365 System needs sufficient space for a computer, printer and selected accessories such as the autosampler. The system dimensions are 491 mm (W) × 538 mm (D) × 258 mm (H) [19.3 in. (W) 21.1 in. (D) 10.1 in. (H)].

WARNING

In view of the electronic components in the Lambda 365 System, it should be positioned well away from any instruments emitting an electromagnetic field and any environment that may cause vibration.

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Power

A grounded AC power supply must be used. The Lambda 365 accepts 100 to 240 VAC at 50 Hz or 60 Hz.

WARNING

In order to protect the system from electrical noise and current surges, please connect the power cable to a GROUNDED OUTLET.

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Environmental Requirements

This system is for indoor use only.

| Requirement | Description |
|--------------------------------|--|
| Indoor use | Altitude up to 2000 m (6,500 ft.) Pollution degree 2 |
| Dimensions | 491 mm (W) × 538 mm (D) × 258 mm (H) [19.3 in. (W) 21.1 in. (D) 10.1 in. (H)] |
| Weight | 23.5 kg (52 lbs) |
| Electrical supply | AC 100-240V, 50/60Hz, 110W |
| Operating Temperature | 15 - 35°C (59 - 95°F) |
| Operating humidity | Less than 80% relative, non-condensing |
| Mains Supply Voltage Functions | ± 10% |

Computer Requirements

| Requirement | Description |
|--------------------|---|
| Hardware | Processor : Intel® Core 1.5 GHz or faster RAM : At least 1GB Hard disk : 50GB with 1GB free Input devices : Mouse and keyboard Monitor : 1024x768 (minimum) Media : CD ROM Drive Port : USB port for the data acquisition |
| Operating System | Microsoft® Windows 7 |
| Output device | Microsoft® Windows compatible printer |

Installation

1. Carefully remove the instrument and all of its components from the shipping container.
2. Verify that all of the components received match the packing list.
3. Verify that all system components required for installation are present, this, must include:
 - Lambda 365 UV-Visible Spectrophotometer
 - Interface cable (USB)
 - Power Cord
 - UV Express Software CD
 - Computer systems (see page 33 for computer requirements)
4. Place the instrument in a location that is compatible with the required environmental conditions for operation. Connect the power cord to the plug; the power plug **must** be grounded.

WARNING

The electric power could be dangerous. Before connecting the power cord, verify that the main power switch is off.

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Unpacking

Heavy Weight

The Lambda 365 weighs 23.5 kg. To avoid injury to personnel or damage to equipment, always use two or more people when lifting or carrying the instrument. NEVER attempt to lift the instrument alone.

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
UV Express Software
Installation


How to Install the UV Express Software

NOTE: *The Sleep Mode in the Power Saver is not recommended to us while running the UV Express software. Set the Power Saver as follows:*


- a. Select **Control Pane> All Control Panel Items > Power Options** from the **Start** menu in the Windows operating system.
- b. Click **Power saver** and click **Change plan settings**.


- c. Click the **Put the computer to sleep** drop-down list and then click **Never** in the Edit Plan Settings window. Click **Save Changes**.


1. Insert the UV Express installation CD into the CD-ROM drive. InstallShield® Wizard starts automatically.
2. InstallShield® Wizard does not start automatically, select **Computer** → **CD Drive**. Double click on **setup.exe**.


3. UV Express InstallShield® Wizard starts. Select **Next**.


4. Enter customer information, if necessary, and select **Next**.


5. The software will be installed in C:\UV Express as a destination folder. To install the software to a different folder, select **Change** and select the desired folder.


6. Verify the current settings and select **Install**.

7. The following window shows the status of the UV Express installation.


8. After the setup is completed successfully, select **Finish**. The UV Express software icon will appear on the Desktop.


Interface Setup

USB Interface Setup


1. Connect the USB cable between the computer and the instrument.
2. Turn on the computer and the instrument.
3. Click **Computer** → **CD Drive** → **Open**.


4. Select **USB Driver** folder.


5. Double click **USBCom-CDM_20828**.


6. Select **Extract**.


7. Select **Next**.

8. The following dialog box will appear. After installation is completed successfully, select **Finish**.


9. Change the USB serial Ports after finishing driver setup. See the following chapter for that procedure.


Setting USB Serial Ports

Procedure


1. Select **Computer** → **Properties**.


2. Click on **Device Manager**.


3. Click on the **Ports (COM & LPT)** to expand the listing. These are the devices currently connected to the COM ports. **USB Serial Port (COMx)** is visible when the driver installation is completed successfully.


4. Double click **USB Serial Port (COMx)** of the Ports (COM & LPT) section.


5. Select on the **Port Settings** tab and click on the **Advanced...** button.


6. Change the parameter values as shown below.


- Verify the parameters and click **OK**.

- Launch the **UV Express** software.
- If the Lambda 365 fails to communicate with the PC, change the COM Port Number as following procedure.
- Open **Advanced Setting for COMx** window following on page 48~50.
- Select the **COM Port number** list to expand it and change the COM port number another one which is not in use from COM 1 to COM 10.

12. Make sure that the changed COM Port Number is applied, and select **OK**.

13. Restart the computer after finishing driver setting.

Operating Procedure


Procedure

- Turn on the power switch and allow the instrument to warm-up for at least 20 minutes.
 - * Check the LED light, which is the lower left hand in front of the system.
 - Power: Blue -> Power on
 - Ready: White-> the communication of between PC and system
 - System: White-> during the checking of system status


CAUTION *Do **not** execute UV Express before finishing the initialization of grating.*

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- Double-click **UV Express** folder.
- Execute **Scan** mode icon in the UV Express Software folder. The following message box will be displayed. Empty the cell holder and close the lid firmly. Click **OK**.

- Start **System Self Test**.

- * **System Self Test** is the procedure that the system checks instrumental conditions by itself.
- * **System Self Test** is performed whenever the system's power is turned on and the software is executed.

CAUTION *Do **not** open the cover during System Self Test.*

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Accessories

Accessories for the Lambda 365

Every day you can count on PerkinElmer to provide you with solutions that deliver reliable performance, control operating costs and maximize operational time. Our complete portfolio of consumables, parts, supplies, training and service helps you meet both routine and demanding measurement challenges. We invest heavily in testing and validating our products to ensure you receive guaranteed compatibility and performance-on-time, every time, for every instrument in your laboratory.

Supplies, replacement parts, and accessories can be ordered directly from PerkinElmer, using the part numbers quoted in the guides provided with the instrument.

See our website:

<http://perkinelmer.com>

PerkinElmer's catalog service offers a full selection of high-quality supplies.

To place an order for supplies and many replacement parts, request a free catalog, or ask for information:

If you are located within the U.S., call toll free 1-800-762-4002, 8 a.m. to 8 p.m. EST. Your order will be shipped promptly, usually within 24 hours.

If you are located outside of the U.S., call your local PerkinElmer sales or service office.

For more detailed information about the accessories, please refer to the instruction sheet for each accessory

Single Cell Holder

Designed for a standard cuvette with a 10 mm pathlength. Optional stirring is available with a separately purchased stir controller.


Water Jacketed Single Cell Holder

Use this accessory to control the temperature of a sample in a standard 10 mm rectangular cuvette. This accessory needs a circulator for control temperature with a separately purchased circulator. Optional stirring is available with a separately purchased stir controller.


8 Position Multicell Holder / Water Jacketed 8-Position Multicell Holder

Use this accessory to measure eight samples automatically. Automatic Water Jacketed 8-Position Multicell Holder accommodates a temperature by circulator. These accessories are controlled by UV Express software.

Micro Cell Holder

Use this accessory to test micro volume samples.

Variable Pathlength Cuvette Holder

Use this accessory to hold rectangular cells with 5, 10, 20, 40, 50, 100 mm pathlengths. A moveable block secures cells in position.


Advanced Transmission Holder

Use this accessory to measure the transmittance spectrum of various solid samples including films, glass, filters, and optical components.

Film Holder

Use this accessory to measure the transmittance spectrum of thin samples such as films and filters.


Variable Angle Transmission Holder

Use for measuring the transmittance of various features of Film, Glass, and Plate, etc according to the variable angle.

50 mm Transmission / Reflectance Sphere

Used for measuring the reflectance or total transmittance of solid samples such as; metals, glass, textiles, paper, gems, and powders.

Fixed Angle Reflectance Holder

Use this accessory to test coatings and semi-conductor wafers; especially effective for measuring specular reflectance.


Single Cell Peltier Holder, Controller

Control the temperature of a single 10 mm pathlength rectangular cell between -5 and 100°C. Peltier control offers fast ramping with recirculating liquid for exceptional temperature control. Integrated purging capability can be used to decrease condensation. This accessory includes a temperature probe, temperature controller, Peltier accessory and all accessories necessary for its operation and connection to the instrument. This accessory is controlled by UV Express software.

with heated reference

without heated reference

Peltier Controller

6-Position Peltier Controlled Cell Changer, Controller

Control the temperature of up to six, 10 mm pathlength rectangular cell between -5 and 100°C. Peltier control offers fast ramping with recirculating liquid for exceptional temperature control. Integrated purging capability can be used to decrease condensation. This accessory includes two temperature probes, a temperature controller, Peltier Temperature Control Unit and all accessories necessary for its operation and connection to the instrument. This accessory is controlled by UV Express software.

with heated reference

without heated reference

Peltier Controller

Nano Stick

Use this accessory to measure micro-volume sample such as DNA, RNA and Protein.

Autosampler

Automate the measurement of up to 360 samples at a time with an autosampler accessory. Two and four rack configurations are available. This accessory system, including the automated sipper, is fully controlled by UV Express software.

Lamp Replacement


D2 Lamp Replacement

1. Make sure that the main power is turned off. If the system and lamp(s) were turned on at least 20 minutes should pass before the lamp changing process is started. This time allows the lamps to cool and return to room temperature.


CAUTION *Use appropriate disposable gloves during the lamp replacement process. Also avoid touching the glass bulb section of the lamp this can leave oil deposits on the bulb which can drastically reduce lamp life.*

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2. Remove the white bolt covers (four each) using screwdriver and remove the M3*12 screws (four each) using a Phillips screwdriver to open the vent cover.


3. Remove the two M3*6 screws on the sub cover using a Phillips screwdriver.

4. Loosen the two screws using a Phillips screwdriver. Note that there is no need to take the screws out from the lamp module.

5. Carefully lift up the D2 lamp module and disconnect the D2 lamp cable connector.

6. Replace the D2 lamp module with a new one and reassemble the parts in reverse order.

Tungsten Lamp Module Replacement

1. Make sure that the main power is turned off. If the system and lamp(s) were turned on at least 20 minutes should pass before the lamp changing process is started. This time allows the lamps to cool and return to room temperature.

CAUTION

Use the Tungsten Lamp Changing Tool (Part No. N4102048) when replacing the Tungsten Lamp. In addition a pair of appropriate disposable gloves should be worn during the lamp replacement process. Please avoid touching the glass bulb section of the lamp this can leave oil deposits on the bulb which can drastically reduce lamp life.

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
Tungsten Lamp Changing Tool (Part No. N4102048)

2. Remove the white bolt covers (four each) using screwdriver and remove the M3*12 screws (four each) using a Phillips screwdriver to open the vent cover.

3. Remove the two M3*6 screws on the sub cover using a Phillips screwdriver.

4. Put the tungsten lamp changing tool on the tungsten lamp module carefully by aligning the notches over the lamp module securing screws.


5. Fasten the cap screw on the tungsten lamp changing tool to fix the tool to the tungsten lamp module.

6. Using a Phillips screwdriver to loosen the two screws that are located on the tungsten lamp module.

7. Carefully lift up the tungsten lamp changing tool which holds the tungsten lamp module and disconnect the tungsten lamp cable connector.

8. Loosen the cap screw on the tungsten lamp changing tool (Part No. N4102048) and separate the tungsten lamp from the tool.

9. Replace the tungsten lamp with a new one and reassemble the parts in reverse order.

Index

| | | | |
|--|----|--|----|
| A | | L | |
| About This Manual | 7 | Laboratory Preparation..... | 32 |
| Accessories | 56 | Lamp Replacement | |
| C | | D2 Lamp..... | 64 |
| Cleaning | 28 | Tunsten Lamp..... | 66 |
| Cleaning the Instrument | 28 | O | |
| Computer Requirements | 33 | Operating Conditions..... | 22 |
| Conventions | | Operating Procedure..... | 48 |
| Note, attention et avertissement | 15 | P | |
| Notes, cautions and warnings | 15 | Precautions | 21 |
| text | 15 | S | |
| D | | Safety | |
| Decontamination | 28 | Summary | 11 |
| E | | Safety Information..... | 21 |
| Electrical Safety | 23 | electrical safety | 23 |
| EMC Compliance..... | 25 | Shipping Kit | 14 |
| Environmental Requirements..... | 33 | U | |
| I | | UV Express Software Installation | 36 |
| Installation | | | |
| Lambda 365..... | 34 | | |
| UV Express Software..... | 36 | | |