

## Thermal Analysis



### Preparation Checklist

- Site Requirements
- Peripherals and Accessories
- Electrical Requirements
- Gas Requirements
- Environmental Requirements
- Safety Requirements
- PC Configuration
- Installation Overview
- Miscellaneous

## DSC 8000: Site Readiness Instructions

### Order Overview

Please review this order. Record any discrepancies between the PerkinElmer order and your Purchase Order, along with any agreements or commitments made by your PerkinElmer Sales Representative that are NOT listed on the order. Let your Customer Care Representative now about these discrepancies and/or commitments with your Site Readiness Confirmation.

### Site Requirements:

#### Bench Space

| Instrument           | Dimensions         |                    |                      |                      |
|----------------------|--------------------|--------------------|----------------------|----------------------|
|                      | Width              | Depth              | Height               | Weight               |
| DSC 8000             | 18 in<br>(45.7 cm) | 25 in<br>(63.5 cm) | 11.5 in<br>(29.2 cm) | 48 lbs.<br>(21.8 kg) |
| DSC 8000<br>w/AS     | 18 in<br>(45.7 cm) | 25 in<br>(63.5 cm) | 19.5 in<br>(49.5 cm) | 48 lbs.<br>(21.8 kg) |
| Computer<br>(approx) | 7 in<br>(17.5 cm)  | 17 in<br>(44.0 cm) | 15.8 in<br>(40.2 cm) | 24 lbs.<br>(11.0 kg) |

A bench top of 100 cm (39.4 in) W x 76.2 cm (30 in) D will accommodate a DSC 8000 system with no accessories.

An autosampler system requires additional clearance above the DSC. 36 in/91.4 cm is required for the autosampler dust cover clearance.

With accessories (Intracooler 2 V/3 V, CLN2 Liquid Nitrogen Cooling System, Printer, Plotter,) additional bench space will be required.

*NOTE: With all cooling accessories a coolant transfer line attaches to the left side of the DSC 8000. Consideration must be given for location of dewar relative to the DSC 8000. Typically, the DSC is located at the left end of the bench. Floor space required for the dewar is approximately 60 cm (24") square. The dewar is 120 cm (48") tall and may NOT fit under the lab bench.*

## Peripherals and Accessories:

| Accessories      | Dimensions         |                     |                      |                  |
|------------------|--------------------|---------------------|----------------------|------------------|
|                  | Width              | Depth               | Height               | Weight           |
| Water Circulator | 8.25 in<br>(21 cm) | 15.75 in<br>(40 cm) | 22.5 in<br>(57.1 cm) | 70 in<br>(32 kg) |
| Intracooler 2    | 15.0 in<br>38.1 cm | 21.0 in<br>53.3 cm  | 29.0 in<br>73.7 cm   | 90 lb<br>40.8 kg |
| CLN2             | 15.0 in<br>38.1 cm | 21.0 in<br>53.3 cm  | 29 in<br>73.7 cm     | 90 lb/<br>0.8 kg |

## Electrical Requirements:

| Power Consumption |                    |
|-------------------|--------------------|
| DSC 8000          | 400 Watts Maximum  |
| Computer          | Not Available      |
| Circulator        | 1440 Watts Maximum |
| Intracooler 2     | 1440 Watts Maximum |
| CLN2              | 150 Watts Maximum  |

| Power Specifications |                                       |
|----------------------|---------------------------------------|
| DSC 8000             | 120 VAC, 2.0 Amps or 240 VAC, 1.1A    |
| Computer             | Not Available                         |
| Circulator           | 120 VAC, 12 Amps or 240 VAC, 7 Amps   |
| Intracooler 2        | 100/115 VAC, 12 Amp or 230 VAC, 6 Amp |
| CLN2                 | 120 VAC, 1.5 Amp or 240 VAC, 1.0 Amp  |

This equipment is designed to operate within 10% of the selected line voltage (Except 240VAC +6%, -10%).

The supply must be smooth, clean and free of transient voltages over 40 volts.

Earth grounding: less than 1 ohm resistance between the grounds of any two components of the system.

| Power Outlets          |                        |
|------------------------|------------------------|
| DSC 8000               | 1 standard outlet      |
| Additional Accessories | 1 separate outlet each |

All outlets should share a common earth ground.

## Gas Requirements:

All gasses and regulators MUST be supplied by customer.

Gas dew point must be lower in temperature than the minimum temperature of the cooling accessory.

If splitting a single gas source is necessary.

Please contact your local service representative for information on the regulator install kit Part number N519-0462. Please visit our website to obtain the address of your nearest PerkinElmer office.

All systems require a "Sample Gas" and "System Gas."

Sample gas flow is controlled by the mass flow controller. This is the gas that is purging the sample cup area. Two inputs are provided, Sample Gas A and Sample Gas B to allow purge gas switching experiments. The flow rate is controlled through the instrument control software.

System gas flow is controlled with a frit (restrictor) system built into the analyzer. The system gas is used to purge the environmental control system of the analyzer. This is the area under and around the sample holder including the area where the cooling accessory is mounted. Flow rates in this area will vary depending on the state of the analyzer. While the cover is in the closed position a lower flow rate is utilized. With the cover open a system boost purge is activated at 8 l/min. The DSC 8000 shipping kit contains a Drier Accessory Kit with replaceable drier cartridge (P/N 0992-0018 Reading Drier Cartridge) and a length of copper tubing. The drier must be installed < 6 feet from the analyzer. If additional copper tubing is required it is the responsibility of the customer to provide

| Gas  | Pressure              | Flow         | Purity Minimum | Dew Point     |
|--|-----------------------|--------------|----------------|---------------|
| Operation at or Above Ambient Temperatures                   |                       |              |                |               |
| <b>Sample Gas:</b><br>Argon, nitrogen, helium, air or oxygen | 20-40 psi/<br>2-3 bar | 20-40 cc/min | 99.95          | <-20 °C       |
| Sub-ambient operation using Intracoolers 2                   |                       |              |                |               |
| <b>Sample Gas:</b><br>Argon, nitrogen, helium, air or oxygen | 20-40 psi/<br>2-3 bar | 20-40 cc/min | 99.95          | <-70/<-100 °C |
| <b>System Gas:</b><br>nitrogen, dry                          | 20-40 psi/<br>2-3 bar | 0.6– 8 l/min | 99.95          | <-70/<-100 °C |
| Sub-ambient operation using CLN2                             |                       |              |                |               |
| <b>Sample Gas:</b><br>Helium, exceptionally dry              | 20-40 psi/2-3 bar     | 20-40 cc/min | 99.95          | <-190 °C      |
| <b>System Gas:</b><br>Nitrogen, dry                          | 20-40 psi/2-3 bar     | 0.6– 8 l/min | 99.95          | <-190 °C      |
| <b>CLN2 Pressure source:</b><br>Nitrogen                     | 20-50 psi/2-3 bar     | N/A          | 99.95          | <-190 °C      |

## Environmental Requirements:

### Laboratory Environment

Clean and dust-free.

Indoor use only on level, vibration-free work surface.

Intracooler must not be located in an enclosed area.

## Safety Requirements:

### Gas Cylinders and Gas Delivery Lines:

Lock down straps should be present on all gas cylinders.

### Ventilation:

Do not operate the Differential Scanning Calorimeter in an enclosed environment without adequate ventilation.

The liquid nitrogen system of PerkinElmer's Differential Scanning Calorimeter emits a small amount of nitrogen during normal operations. If your laboratory is not properly ventilated, the level of oxygen in the laboratory may fall below the normal range. Please follow the applicable laboratory ventilation standards to ensure that an appropriate oxygen level is maintained.

### PC Configuration:

Pyris Software is localized to only English and Japanese. The user must set the locale to either of these two settings.

Due to numerous differences in PC hardware, PerkinElmer cannot guarantee that our software will run on a customer-supplied computer.

PerkinElmer installation of a customer-supplied computer is available for an additional fee.

PerkinElmer is not responsible for problems caused by unspecified system components, software, and/or accessories. A maximum of one hour is allowed for installation of a computer and software of a non-PerkinElmer supplied computer. The additional time it takes to verify this type of problem is billable at the current service rate. It is advisable the customer's IT support be available as needed.

### Software:

As per requirement.

### Installation Overview:

On receipt of the equipment, the customer is to inspect the packaging for physical damage. If damage is present the shipping container should be opened to verify no physical damage to the instrument has occurred and the customer must notify the shipper immediately.

*NOTE: Unpacking will only be performed by approved PerkinElmer personnel and contents inventoried.*

### Physical Installation: (Instrument Only)

The physical installation will vary based on system configuration.

### Physical Installation: (Accessories)

As required.

### Installation Test Standards:

Our Service Engineer will test the instrument in order to insure that it meets the functionality requirements.

## Miscellaneous:

### Cooling System Requirements:

In order to get the maximum instrument performance, it is recommended that a cooling system be used for DSC operations at or above ambient temperature.

### For Operation at or Above Ambient Temperature:

The DSC 8000 can be configured with a turbulent chamber.

Using the turbulent chamber the DSC 8000 can be configured with a tap water circulating system or circulating system (PolyScience Chiller).

With the tap water configuration (with drain) or a circulating device (chiller), either must be located conveniently near the DSC 8000.

### For Subambient Operations:

It is necessary to install an Intracooler or a CLN2 accessory.

### Sample Preparation:

All reference materials required for installation check-out are shipped in the instrument start-up kit.

### If IQ/OQ Validation is Required for This Instrument:

Please contact your Customer Care Representative or visit our website to obtain the address of your nearest PerkinElmer office.

## Miscellaneous Accessories and Spares

### Upgrade a DSC 8000 to a DSC 8500

| Part Number | Description                        |
|-------------|------------------------------------|
| N5340509    | Upgrade kit for DSC8000 to DSC8500 |
| N5320129    | Balance Tare Weight, Medium        |

### Crimper Presses for DSC 8000/8500

| Part Number | Description   |
|-------------|---|
| Part Number | Description   |
| 02190048    | Standard Sample Pan Crimper Press (for standard pans) |
| B0139005    | Universal Crimper Press (for robotic pans)            |
| 03190030    | Sample Holder Cover Reforming Tool                    |

## Pans and Covers for DSC 8000/8500

### Aluminum Open/Crimped Pans

| Part Number | Description  |
|-------------|--|
| 02190041    | Standard Aluminum Pans & Covers (400/Pkg)                            |
| B0198030    | Aluminium Pans and quartz disk Covers for Photocalorimeter (100/Pkg) |
| B0196858    | Aluminium Pans for Photocalorimeter (100/Pkg)                        |
| B0181091    | Quartz disks for Photocalorimeter pans (100/Pkg)                     |

### Aluminum Vented/Pierced Pans

| Part Number | Description   |
|-------------|---|
| B7001014    | Aluminum covers – pierced (400/Pkg)                                   |
| N5190788    | Volatile Aluminum Pans and Covers, 20 micro liter – pierced (100/Pkg) |
| B0143018    | Aluminum (vented), 30 micro liter (400/Pkg)                           |
| B0143019    | Aluminum (vented), 50 micro liter (400/Pkg)                           |

### Aluminum Hermetically Sealed/Volatile Pans

| Part Number | Description  |
|-------------|--|
| B0143015    | Aluminum Pans, 10 micro liter (400/Pkg)            |
| B0143016    | Aluminum Pans, 30 micro liter (400/Pkg)            |
| B0143017    | Aluminum Pans, 50 micro liter (400/Pkg)            |
| 02190062    | Aluminum Pans and Covers, 20 micro liter (400/Pkg) |
| B0169319    | Aluminum Pans and Covers, 10 micro liter (400/Pkg) |
| B0169320    | Aluminum Pans and Covers, 30 micro liter (400/Pkg) |
| B0169321    | Aluminum Pans and Covers, 50 micro liter (400/Pkg) |

### High Pressure Sample Pans

| Part Number | Description   |
|-------------|---|
| 03190218    | Large Volume Stainless Steel Pans, Covers and O-rings, 60 micro liter, 24 Atm (20/Pkg)            |
| 03190029    | Large Volume Stainless Steel Pans, Covers and O-rings, 60 micro liter, 24 Atm (1000/Pkg)          |
| B0182901    | Re-useable Stainless Steel Pans, Covers and 20 gold-plated seals, 30 micro liter, 150 Atm (5/Pkg) |
| B0182902    | Gold Pans and Covers, 30 micro liter, 150 Atm (5/Pkg)   |
| B0182903    | Titanium Pans and Covers, 30 micro liter, 150 Atm (5/Pkg)   |

### Aluminum sample pan kit

| Part Number | Description                                    |
|-------------|--|
| B0510800    | Aluminum Pans and Covers Starter Kit (400/Pkg) |

### Specialty Sample Pans and Covers

| Part Number | Description                                   |
|-------------|---|
| N5203115    | HyperDSC™ Aluminum Sample Pans (100/Pkg)      |
| N5190180    | Alumina Sample Pans and Covers (6/Pkg)        |
| 02190042    | Standard Gold Sample Pans and Covers (10/Pkg) |
| 03190025    | Graphite Sample Pans and Covers (4/Pkg)       |

## Furnace Covers for DSC 8000/8500

| Part Number | Description  |
|-------------|--|
| 03190031    | Vented Platinum Furnace Covers with 2 holes (6/Pkg)    |
| 04190299    | Vented Platinum Furnace Covers with 2 holes (2/Pkg)    |
| 03190032    | Unvented Platinum Furnace Covers with no holes (6/Pkg) |
| B0139122    | Autosampler Furnace Covers (2/Pkg)                     |
| B0182981    | Large Autosampler Furnace Covers for Tall Pans (2/Pkg) |

Refer to the consumables and accessories reference catalogue for full descriptions of pans and covers.

## Calibration Supplies

| Part Number | Description  |
|-------------|--|
| 03190033    | Indium calibration reference   |
| 03190034    | Tin calibration reference  |
| 03190035    | Lead calibration reference   |
| 03190036    | Zinc calibration reference   |
| 02190045    | Encapsulated calibration reference kit with large and small indium, tin and lead |
| N5190762    | Encapsulated calibration reference kit with indium and zinc                      |

PerkinElmer, Inc.  
940 Winter Street  
Waltham, MA 02451 USA  
P: (800) 762-4000 or  
(+1) 203-925-4602  
[www.perkinelmer.com](http://www.perkinelmer.com)



For a complete listing of our global offices, visit [www.perkinelmer.com/ContactUs](http://www.perkinelmer.com/ContactUs)

Copyright ©2016, PerkinElmer, Inc. All rights reserved. PerkinElmer® is a registered trademark of PerkinElmer, Inc. All other trademarks are the property of their respective owners.