# APPLICATION NOTE



# Gas Chromatography/ Mass Spectrometry

Author: Sharanya Reddy PerkinElmer, Inc. Shelton, CT

Rapid Fingerprinting of Volatile Components in Coffee Varieties Using the Portable Torion T-9 Portable GC/MS System

## Introduction

Coffee is widely consumed as a beverage because of the stimulating effect it produces in humans. The aroma of coffee contributes to the flavor and taste of the

beverage and has consequently led to extensive research on coffee. Hundreds of volatile/semivolatile aroma and flavor compounds have been identified in coffee using traditional laboratory based headspace GC/MS systems<sup>1</sup>. We present a study for rapid fingerprinting of coffee volatile/semi-volatile compounds using solid phase micro-extraction (SPME) coupled to a portable GC/MS system for separation and detection.

The on-site analysis of coffee using portable technology can be used for quick quality control check of raw and finished products, comparison of competitor products, analysis of storage conditions or for process development.



# Experimental

## **Sample Preparation**

Coffee (1 gm) was placed in headspace vials (20 mL vial volume), capped and placed at room temperature for at least an hour to allow for saturation of the aroma volatiles in the headspace. The SPME fiber was directly exposed to the headspace vapors in the vial for 15 seconds prior to analysis.

#### Table 1. GC and MS conditions.

| Sampling             | SPME   |
|----------------------|--|
| SPME Phase           | Divinylbenzene/Polydimethylsiloxane<br>(DVB/PDMS, 65 μm) |
| GC Injector Temp     | 270°C  |
| GC Column:           | MTX-5, 5 m X 0.1 mm, 0.4 µ df                            |
| GC Carrier Gas       | Helium, 0.2 ml/min                                       |
| GC Column Temp       | 50-270 °C at 2 °C/sec, end hold time for 60 sec          |
| Transfer Line        | 250 °C   |
| Injector Split ratio | 10 to 1  |
| Mass analyzer        | Toroidal ion trap  |
| Mass Range           | 42-500 Da  |
| Detector             | Electron multiplier                                      |

#### **Results and Discussion**

Triplicate analysis of coffee (variety #1) is shown in Figure 1. The overlay of the Total Ion Current (TIC) suggests the analysis is very reproducible between injections.



*Figure 1.* Overlay of the TIC of triplicates analysis of headspace vapors of coffee sample (variety #1).

The decaffeinated form of variety #1 coffee was similarly analyzed in triplicates (Figure 2) and also showed excellent reproducibility of analysis.



Figure 2. Overlay of the TIC of triplicates analysis of headspace vapors of decaffeinated coffee sample (variety #1).

The overlay of TIC of caffeinated and decaffeinated samples of coffee (variety #1) showed peaks with similar retention time but in many cases with varying intensities suggesting similar compounds are present in the two samples but at different concentrations (Figure 3). An extra peak was observed in the caffeinated sample that was not observed in the decaffeinated coffee.



Figure 3. Overlay of the TIC of caffeinated and decaffeinated coffee variety #1 samples.

The extraneous peak observed in the caffeinated coffee resembled the spectra of toluene and comparisons confirmed its identification (Figure 4). The presence of toluene in the caffeinated sample is not surprising as literature suggests toluene can be produced in roasted coffee<sup>2</sup>.



*Figure 4.* Mass spectra of peak observed at retention time at 32.3 sec in coffee sample (variety #1, Figure 3) matched with spectra of toluene (head to tail comparison) with a match factor of 98.5.

Two other varieties of coffee (variety #2 and #3) were analyzed and the TIC overlay of three coffee types (Figure 5) showed some distinct profiles between the samples. For instance, a peak at ~30 secs was common to #1 and #2 and not observed in #3. The spectra of this peak, was matched against NIST library and the probability of the match suggested it was likely methyl pyrrole.



Figure 5. Overlay of the TICs of three different coffee varieties.

Literature studies suggest that methyl pyrrole degradation is accelerated in coffee with high moisture content<sup>3</sup>, which may be the reason for variety #3 having little or no methyl pyrrole.

### Conclusion

A brief study of coffee fingerprinting using the Torion T-9 coupled to SPME sampling resulted in reproducible profiles for triplicate analysis. Separation of the volatile/semi-volatile compounds by GC and the specificity of detection provided by MS along with spectral matching to the NIST library helped identify differences in coffee varieties. The Torion T-9 portable GC/MS provides a rapid and reliable screening technology to study coffee aroma profiles within a matter of a couple of minutes making it an ideal tool for on-site analysis.

#### References

- Volatile compounds in Food. Nijssen L. M., Visscher C.a., Maarse H., Willemsense L. C., Boelens, M.H., 7th ed.; TNO Nutrition and Food Research institute: Zeist, The Netherlands.
- 2. Cancer prevention and the environmental chemical distraction. Ames B., Gold L.S., Politicizing Science: the Alchemy of Policymaking. Micheal Gough ed., Hoover Institute Press, Stanford, California. pp 117-142.
- Quenching method, moisture content and aroma stability of roast and ground coffee. Baggenstoss J., Escher F, Expression of Multidisciplinary Flavor science – 12th Weurman Symposium', Zurich University of Applied Science, Wintherthur, 2010, Chapter 6, pp 271-274.

PerkinElmer, Inc. 940 Winter Street

940 Winter Street Waltham, MA 02451 USA P: (800) 762-4000 or (+1) 203-925-4602 www.perkinelmer.com



For a complete listing of our global offices, visit www.perkinelmer.com/ContactUs

Copyright ©2016-2017, PerkinElmer, Inc. All rights reserved. PerkinElmer® is a registered trademark of PerkinElmer, Inc. All other trademarks are the property of their respective owners.

PKI