
Figure 1 The FIAS-400 Flow Injection System.

The On-line Preconcentration Kit for the FIAS-400 (or 200) Flow Injection System provides automated analyte preconcentration capabilities which significantly extend the fields of application for flame atomic absorption.

Some of the key benefits you can obtain with the On-line Preconcentration kit for FIAS include:

- Analyte preconcentration for improved detection limits in your sample matrices.
- Complete, computer-controlled automation for operational simplicity.
- Significantly reduced sample preparation time.
- Reduced risk of contamination during sample preparation.
- Improved accuracy by separating potentially interfering matrix components prior to the analytical measurement.


Figure 2 FI Preconcentration Tubing Installation using sorbent extraction.

With many sample types, the analyte concentration can be too low for direct determination by flame atomic absorption. Frequently, time-consuming sample preparation techniques are required to extract or preconcentrate the elements of interest. With the On-line Preconcentration Kit for FIAS such procedures can now be fully automated using advanced flow injection techniques. On-line preconcentration of trace elements enhances the detection limits of flame atomic absorption. For many applications, detection limits in your samples can be improved significantly.

Two different preconcentration techniques are currently offered. The most commonly used method is on-line preconcentration with sorbent extraction. The sample is mixed automatically with a chelating agent and passed to a specially designed column containing C₁₈-coated silica or other sorbent materials. After the desired preconcentration time, adsorbed metal chelates are eluted with alcohol and transported to the flame via the nebulizer. All steps, including sampling, mixing, preconcentration on the column and elution, are fully automated and controlled by the AA instrument's computer.


Figure 3 Determination of 3 µg/L of Cd in a solution containing 10 % dissolved solids after 40 seconds preconcentration.

With proper selection of a suitable chelating agent, column material and eluent, this technique provides very selective preconcentration of specific metals such as Cu, Cd, and Pb. This allows them to be separated from undesirable matrix components, such as high concentrations of alkali metal or alkaline earth salts. For example, using automated on-line preconcentration, these elements have been successfully determined directly in drinking water and coastal seawater using flame atomic absorption.

Ion exchange is another proven method for the FIAS automated preconcentration system. For example, this technique has been used for the determination of gold in ores using an XAD-8 column. A preconcentration factor of 15-20 can easily be obtained within a preconcentration time of only 20 seconds. This improves the flame AA detection limit for gold in ore to approximately 2.0 µg/L.

Automatic on-line preconcentration saves sample preparation and analysis time and costs. It also enhances sample throughput and reduces the potential for interferences for better analytical accuracy.

Determination of trace elements in NBS 1643b water sample

Element	Measured Value (µg/L)	Certified Value (µg/L)
Cu	22 ± 0.4	22.3 ± 0.4
Pb	23 ± 0.3**	23.7 ± 0.7
Cd	19 ± 0.3	20 ± 1

Figure 4 Determination of Cu, Pb, and Cd in a certified water sample.

Ordering Information

B050-8593

Preconcentration Kit for FI-Flame AA using the FIAS-400. For automatic on-line preconcentration with flame AA. Includes selected tubing and connectors. Requires but does not include appropriate preconcentration columns (see below).

B050-1593

Preconcentration Kit for FI-Flame AA using the FIAS-200. For automatic on-line preconcentration with flame AA. Includes one 5-port sampling/switching valve, tubing and connectors. Requires but does not include appropriate preconcentration columns (see below).

B050-4047

C₁₈ Preconcentration Column. PTFE column packed with C₁₈ material for use with preconcentration kit B050-1593. Pack of 10.

B050-4048

CAD-8 Preconcentration Column. PTFE column packed with XAD-8 material for use with preconcentration kit B050-1593. Packet of 10.

ALL ANALYTICAL INSTRUMENTS AND SYSTEMS
MANUFACTURED BY PERKIN-ELMER ARE CERTIFIED
UNDER THE QUALITY REQUIREMENTS OF ISO 9001.

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Perkin-Elmer Corp., Analytical Instruments, Main Ave. (MS-12), Norwalk, CT 06856 U.S.A. Tel: (203) 762-1000
Bodenseewerk Perkin-Elmer GmbH, Postfach 101164, W-7770 Überlingen, Bundesrepublik Deutschland, Tel: (07551) 81-0
Perkin-Elmer Ltd., Post Office Lane, Beaconsfield, Bucks HP9 1 QA, England Tel: Beaconsfield (0494) 676161