

Explosives by HPLC Analysis per U.S. EPA Method 8330

Author:
Vivek Dhyani

PerkinElmer Life and
Analytical Sciences
710 Bridgeport Avenue
Shelton, CT 06484

Introduction

U.S. Environmental Protection Agency (EPA) Method 8330 is intended for the trace analysis of explosive residues by high performance liquid chromatography (HPLC) using a UV detector. This method is used to determine the concentration of the following compounds in a water, soil, or sediment matrix:

- | | |
|----------------------------|-------------------------------|
| 1. HMX | 8. 2-Amino-4,6-Dinitrotoluene |
| 2. RDX | 9. 4-Amino-2,6-Dinitrotoluene |
| 3. 1,3,5-Trinitrobenzene | 10. 2,6-Dinitrotoluene |
| 4. 1,3-Dinitrobenzene | 11. 2,4-Dinitrotoluene |
| 5. Tetryl | 12. 2-Nitrotoluene |
| 6. Nitrobenzene | 13. 4-Nitrotoluene |
| 7. 2, 4, 6-Trinitrotoluene | 14. 3-Nitrotoluene |

Experimental

The analysis was carried out using a PerkinElmer® Series 200 HPLC System, including Binary Pump, Autosampler and Photo Diode Array (PDA) Detector. The column used was the Brownlee™ Validated C18 column, 250 X 4.6 mm, 5 µm. The mobile phase conditions were 50:50 methanol/water, with a flow rate of 1.5 mL/min, run at ambient temperature. The UV/Vis detector was set to 220 nm and the injection volume was 10 µL. The standard solution was diluted in acetonitrile before injection.

Results

The chromatographic separation is shown in Figure 1. The peaks are numbered corresponding to the order provided in the introduction, above.




Figure 1: Chromatographic separation of 14 explosives by reversed-phase chromatography. The HPLC method conditions and component identification are provided above.

Conclusion

Using the RP C18 column, baseline separation of 8 explosives is possible. In addition, 2-Amino-4,6-Dinitrotoluene and 4-Amino-2,6-Dinitrotoluene are resolved, though not baseline separated. Tetryl coelutes with nitrobenzene and 2,6-dinitrotoluene coelutes with 2,4-dinitrotoluene.

PerkinElmer, Inc.
940 Winter Street
Waltham, MA 02451 USA
Phone: (800) 762-4000 or
(+1) 203-925-4602
www.perkinelmer.com



For a complete listing of our global offices, visit www.perkinelmer.com/lasoffices

©2007 PerkinElmer, Inc. All rights reserved. The PerkinElmer logo and design are registered trademarks of PerkinElmer, Inc. Brownlee is a trademark and PerkinElmer is a registered trademark of PerkinElmer, Inc. or its subsidiaries, in the United States and other countries. All other trademarks not owned by PerkinElmer, Inc. or its subsidiaries that are depicted herein are the property of their respective owners. PerkinElmer reserves the right to change this document at any time without notice and disclaims liability for editorial, pictorial or typographical errors.

The data presented in this Field Application Report are not guaranteed. Actual performance and results are dependent upon the exact methodology used and laboratory conditions. This data should only be used to demonstrate the applicability of an instrument for a particular analysis and is not intended to serve as a guarantee of performance.