

## ***FL 6500/8500 Fiber Optic Probe Installation Instructions***

This instruction sheet describes the installation of this accessory which is used with the FL 6500/8500 Fluorescence Spectrometer.

**NOTE:** *Read these instructions before you install this accessory.*

### ***Contacting PerkinElmer***

Supplies, replacement parts, and accessories can be ordered directly from PerkinElmer, using the part numbers.

See our website:

<http://perkinelmer.com>

PerkinElmer's catalog service offers a full selection of high-quality supplies.


To place an order for supplies and many replacement parts, request a free catalog, or ask for information:

If you are located within the U.S., call toll free 1-800-762-4000, 8 a.m. to 8 p.m. EST. Your order will be shipped promptly, usually within 24 hours.

If you are located outside of the U.S., call your local PerkinElmer sales or service office.

### ***Features***

- Sample measurement with simultaneous Excitation and Emission beam reflections
- Various sample measurement for solid, liquid, powder sample, etc.
- Include Probe Holding Stand



**Figure 1 FL 6500/8500 Fiber Optic Probe [P/N:N4201021]**



PerkinElmer, 710 Bridgeport Avenue,  
Shelton, CT 06484-4794, U.S.A

**Produced in the USA.**

## *Dimensions and Specifications*

Physical characteristic		Specification
Fiber Optic Probe Accessory	Fiber length (mm)	2000
	Fiber diameter ( $\mu\text{m}$ )	600
	Probe diameter (mm)	3.14
	Probe length (mm)	100
	Weight (Kg)	0.16
Base Plate	Dimensions (mm)	130 x 267 x 97 (WDH)
	Weight (Kg)	0.72
Probe Holding Stand	Dimensions (mm)	153 x 153 x 145 (WDH)
	Weight (Kg)	1.04

## *Measurable Sample*

Description
a. Liquid sample (Minimum 7 mm inlet diameter sample container, Maximum 70 mm long immersing probe length)
b. Solid sample (Power, Film, Paper, Plastic, No limit to the thickness and diameter)


## *Configuration of the Fiber Optic Probe*



**Figure 2 Fluorescence Probe**




**Figure 3 Stand Module**



**Figure 4 Probe Holding Stand**


## *Installation*

1. Prepare the FL 6500/8500 Fluorescence Spectrometer to install this accessory.
2. Connect the power cord and the communication cable.
3. Loosen the accessory fixing bolt to take apart the existing cell holder.



**Figure 5 Loosening the Accessory Fixing Bolt**

4. Pull out the cell holder by hand.



**Figure 6 Pulling out the Cell Holder**

5. Remove the front cover plate.




**Figure 7 Front Cover Plate**

6. After checking the pogo pin position of the sample compartment, attach the Fiber Optic Probe Base Plate to the pogo pin.



**Figure 8 Install the Base Plate**

7. Tighten the accessory fixing bolt.



**Figure 9 Tightening the Accessory Fixing Bolt**

8. Connect the fiber for excitation and emission.


Excitation Adapter: SPECTROMETER line

Emission Adapter: LIGHT SOURCE line



**Figure 10 Connecting the Fiber**

9. Pass the sample line of Optic Fiber outward from inside the sample compartment.



**Figure 11 Connecting the Fiber**

**NOTE:** *End side of the Fiber Optic is fixed to the holder in order to avoid scratching because they are polished.*

10. Fix the front cover plate attached to the Fiber Optic Probe and tighten the bolts.




**Figure 12 Fixing the Front Cover Plate**

11. Select one measurement method among the following methods depending sample types.

a. Solid Samples

- i) Insert the Probe into the stand module hole and fasten the knob.




- ii) Measure the sample as shown below.



b. Liquid Samples


- i) Insert the Probe into the probe hole of the Probe Holding Stand and fasten the knob. And then, measure liquid and powder samples as shown below.



**NOTE:** When liquid sample is measured, the bottle container should be light blocked.

## Measurement

1. Double click on the **Spectrum FL** software.
2. Check the recognition of Accessory.



3. Click measurement mode.
4. Set up the measurement parameters.

**NOTE:** For more detail of method, refer to *Spectrum FL Software Users Guide*.


5. Click **Save** to save the method after setting up the parameters.
6. Mount the probe according to the sample type and click **Run** icon.
7. Input the sample name and click **OK**.
8. The results are shown the result window.
9. Save or print the data as required.

## Troubleshooting

### *When fluorescence intensity is too low*

1. Check whether the fibers are connected to the right positions, Excitation and Emission.
  - **Distinction of fibers:** The fiber for Excitation has one light path line on the connector tip and the fiber for Emission has a bundle of light path lines.

**Excitation**



**Emission**



### *When Fluorescence is not Measured Through the Fiber Optic Probe*

1. Check whether the fibers are connected to the collimator lens tightly, If they are loose, tighten the connections again.

***When the Reproducibility of Fluorescence Spectrum is Reduced***

1. Holding the probe by hand for measurement is difficult to maintain the stable fixation of the probe tip to the sample. Use the Probe Holding Stand for secured fixation.